Dentsply Denture Offering					
			Design	Guide Codes	
Code Name	Description	Option Name	Description	Option Name	Description
1. Repair Model	Indicates if holes and/or bubbles are going to be repaired	1.1. Dont Repair Model	The model is kept as it is sent Minor holes and bubbles will be filled and	-	
2. Post Dam	Indicates if the case needs post	1.2. Repair Holes and Bubbles	smoothed. (default) Denture base will have post dam		
3. Retention	dam or not Indicates how much block out should be applied to each case	2.2 No post dam	Denture base will not have post dam Indicates to block out any the heavy retention		
		3.1 Minimal Block Out 3.2 Full Block Out	areas Indicates to block out all areas of retention.	-	
		3.3 Anatomical Block Out	Indicated in immediate dentures mostly. Indicates to block out incisive papilla, mid palatine suture and rugae if needed.		
4. Teeth Libraries S. Bite Classes Reference	Indicates which library, shape and posterior anatomy you want to use in each case Indicates what bite class is going to be each case. Using Reference, bite block or immediate	4.1 Denstply Libraries	Indicates which denstply you are going to use	4.1.1. Denstply HC	High Characterized Portrait libraries High Characterized Genios libraries (default) Standard Portrait libraries
		5.1 Reference	Cases with Pre Op Scans (Previous Dentures)	4.1.2. Denstply Standard 5.1.1 Copy Reference Exactly	Standard Genios libraries Teeth size, Teeth positioning is going to be followed exactly.
		5.2 Bite Block	Cases with wax rim blocks (Wax Bite Rims)	5.1.2 Existing vs Projected Bite Class 5.2.1 Follow Bite Block 5.2.2 Follow Patient Ridges	Indicate which bite class is going to be achieved over the patients ridge position Indicate to follow facial position and markings in bite block Teeth arches to be centered in the patient's ridge
		5.3 Immediate	Cases with existing teeth to be removed	5.2.3 Existing vs Proiected Bite Class 5.3.1 Follow existing u/l teeth 5.3.2 Follow existing upper teeth 5.3.3 Follow existing teeth	Indicate which bite class is going to be achieved over the patients ridge position Indicate's that we'll be following upper and lower position as it was in patient mouth. (For dual cases) Indicate's that we'll be following only upper position text, his researding lower. (For dual cases) Indicate's that we'll be following the existing text hin the arch. (For single cases)
6. Posterior Bite Class	Indicates what bite class is going to be selected in the posterior	6.1 Antomical Normal Occlusion 6.2 Lingualized Occlusion 6.3 Crossbite	Cusp-fossae relationship Upper buccal out of occlusion Cusp-fossae relationship inverted	5.3.4 Existing vs Projected Bite Class	Indicate which bite class is going to be achieved over the patients ridge position
7. Teeth Set Up Changes	Indicates changes in teeth positioning with reference to the reference denture OR bite block	7.1 Vertical Dimension (VDO)	Indicates what changes in the patient's current vertical dimension are going to be designed	7.1.1 Maintain VDO 7.1.2 Open VDO 7.1.3 Close VDO	The current vertical dimension is going to be kept as customer send it. The current vertical dimension is going to be open XX milimeters requested by the customer, always from the lower jaw in the PIN section. "Exception: If customer request to open upper jaw OR vertical falls under the customer's responsibility. The current vertical dimension is going to be closed XX milimeters requested by the customer, always from the lower jaw in the PIN section. "Exception: If customer request to close upper jaw OR vertical falls under the customer's responsibility.
		7.2 Lip Support	Indicates if designers follow OR change existing buccal position in reference denture, bite block or existing teeth (Immediates)	7.3.1 Maintain Lip Support	The current lip support position is going to be maintained based on the reference of the patient's existing teeth.
				7.3.2 Outward	The current lip support position in the reference denture, bite block or patient's existing teeth is going to be moved outward by XX milimeters requested by the customer
				7.3.3 Inward	The current lip support position in the reference denture, bite block or patient's existing teeth is going to be moved inward by XX milimeters requested by the customer
		7.3 Occlusal Plane	Indicates if designers follow OR change existing occlusal plane in reference denture, bite block or existing teeth (Immediates)	7.4.1 Maintain Occlusal Plane	The current occlusal plane is going to be maintained based on the reference of the patient's existing teeth.
				7.4.2 Raise Occlusal Plane 7.4.3 Lower Occlusal Plane	The current occusal plane in the reference denture, bite block or patient's existing teeth is going to be raised by XX milimeters requested by the customer The current occusal plane in the reference denture, bite block or patient's existing teeth is going to be lowered by XX
		7.4 Midline	Indicates if designers follow OR change existing midline position in reference denture, bite block or existing teeth (Immediates)	7.5.1 Maintain midline	milimeters requested by the customer The current midline is going to be maintained based on the reference of the patient's existing teeth.
				7.5.2 Move Right	The current midline position in the reference denture, bite block or patient's existing teeth is going to be moved to the right by XX millimeters requested by the customer.
				7.5.3 Move Left	The current midline position in the reference denture, bite block or patient's existing teeth is going to be moved to the left by XX millimeters requested by the customer.
8. Festooning	Indicates what type of gingiva characterization each case is going to have.	8.1 Festooning Type	Indicates how deep the festoon characterization is going to have each case.	8.1.1 None 8.1.2 Soft 8.1.3 Medium	There is no festoons over the teeth. Gingival margins and papillas to be designed only. Flat buccal surface. There is soft festoons over the teeth. There is medium festoons over the teeth.
		8.2 Stippling	Indicates if extra characterization is needed (Buccal)	8.1.4 Heavv	There is heavy festoons over the teeth. Plus frenum/mobile tissue characterization.
9. T-Bars Denture Base	Indicates if extra support is needed in the denture base.	9.1 Two u/I T-bars	Indicates to use the regular print support	]	
		9.2 Two upper bars and three lower bars	Indicates to put a third bar for support in the lower denture base, print with ASIGA.		

**3shape<sup>▶</sup>**Design Services