

# FULLCONTOUR

by 3shape 


We can design the following indications in Exocad as long as the order is set up correctly:

- Crowns
- Copings
- Bridges
- Frame Bridges
- Models
- Veneers
- Inlays/Onlays
- Abutments (Need the exact implant kit put in the instructions box)
- LUCID
- All on X (Limited)

## How to upload an Exocad case to the FullContour platform

### Step1:

Set Up your order and import your scans. Once complete click save.




### Step2:

Click on the 'Open in explorer' button. This will open your current case in the CAD-DATA folder




### Step3:

Click on **CAD-DATA** in the folder path and it will highlight the case you are working on. This folder needs to be zipped up. Nothing should be removed from the folder.




### Step4:

Zip the folder by right clicking, then click on 'Send to' then click on 'Compressed zipped folder'



### Step5:

Open [portal.fullcontour.com](https://portal.fullcontour.com) in your web browser. Drag and drop your compressed zip folder onto the orders tab




**Step 6:**

Fill in the order form to match **EXACTLY** what is selected in the Exocad order.

If you need another line click the + icon on the right

\* **Guide** is where the visual instructions to the designers sit. Always select **Guide** instructions instead of free typing in the instruction field.



The **Destination** area is where FullContour will send the case to be manufactured. If you have a manufacturer turned on on your account, Select the proper one pertaining to the case. In-House Manufacturing means you will download the case and either manufacture it yourself or send it to a manufacturer on your end.

\* Argen manufacturing cannot be chosen for Exocad orders. You will have to select In House Manufacturing and send to the Argen link once the case is completed

**Step 7:**

Once the order is filled out correctly click 'Send New Orders'. Once the order is sent you will see a green confirmation box



**If you need any further assistance please contact customer support**

Email: [support@fullcontour.com](mailto:support@fullcontour.com)

**US Support: 602-688-4133**  
Monday - Friday: 7:00am - 4:00pm (MST)



English Spanish

**International Support: +44 1595 728586**  
Monday - Friday: 9:00am - 5:00pm (CEST)



English Spanish French Russian